
An Addition to the growing "New Generation" of minimal invasive Retinal Implants

- Used to create additional scleral indentation when placed directly over a break.
- Especially useful for small horseshoe shaped tears that "Fishmouth".
- May be used in conjunction with a vitrectomy when broad support is needed and removal of all traction is not desired.
- Fewer components needed on selected cases.
- Designed to enhance the buckling effect of a sleeve when used in a simple encircling procedure.
- Designed for use with Labtician Circling Bands and Labtician Silicone Strips.
- See reverse for surgical details.

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The Wong Sleeve was developed in cooperation with Dr. David T. Wong, Assistant Professor
University of Toronto faculty of Medicine, Toronto, Ontario, Canada

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Using the Buckling Sleeve as a Buckling Element

by David T. Wong, MD, FRCSC

Introduction



The development of the silicone sleeve has made joining solid silicone bands easier. When using a tire element in conjunction with an encircling band, the silicone sleeve is typically placed 180 degrees from the tire element to reduce bulk and conjunctival erosion.

With the development of minimal buckling techniques where the band becomes the buckling element, the sleeve itself can be used to create additional scleral indentation to further support retinal breaks and reduce vitreous traction.

Technique

A scleral belt loop (or suture loop) is created adjacent to the marked retinal break. The sleeve is placed over the marked break to join the 2 ends of the band. After the band is tightened and the desired indentation is achieved, the sleeve can be further adjusted into position. The effect is similar to a low indentation meridional implant.

Further indentation under the sleeve may be created by placing the belt loop over the retinal break. Position the sleeve within the belt loop, in a manner similar to an intrascleral implant. The belt loop must be sufficiently wide and at least 50% of the scleral thickness to distribute the increased forces.

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NOTE: This Poster was in part originally presented at the ASRS 2004 and is reprinted with permission.

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