



MIRA *Ophthalmic Cryo*

An ISO 13485 Company



*The Ophthalmic Cryo
that provides instantaneous
freezing and defrosting
with tip-freeze-only
probes*





Cryo Systems

CR4001 Ophthalmic Cryo System Includes:

- CR4000 Ophthalmic cryo console with footswitch
- CR4160 25 foot exhaust hose
- CR4260 60 inch pressure hose
- CR5350 Sintered metal filters (pack of three)
- CR43xx Gas connector/filter (specify type of gas)

CR4002 Ophthalmic Office Cryo System Includes:

- CR4000 Ophthalmic cryo console with footswitch
- CR4010 (2) 2.5mm curved retinal probe
- CR4160 25 foot exhaust hose
- CR4260 60 inch pressure hose
- CR5350 Sintered metal filters (pack of three)
- CR43xx Gas connector/filter (specify type of gas)

CR4003 Ophthalmic O.R. Cryo System Includes:

- CR4000 Ophthalmic cryo console with footswitch
- CR4010 (2) 2.5 mm curved retinal probe
- CR4150 Dual cylinder cryo carrier
- CR4160 25 foot exhaust hose
- CR5350 Sintered metal filters (pack of three)
- CR42xx (2) Gas connector (specify type of gas)

Cryo Specifications

- Console Dimensions:** Width 5 inches (13.3 cm)
Height 6.5 inches (16.5 cm)
Depth 11 inches (27.9 cm)
- Weight:** 16 pounds (7.3 kg)
- Front Panel Gauge:** Indicates incoming cylinder gas pressure
- Temperature Selector:** -25°C, -40°C, -55°C, tolerance $\pm 5^\circ\text{C}$
- On/Off Switch:** Turns console on/off (controls gas lines)
- Footswitch:** Controls freezing operation
(depress to freeze, release to defrost)
- Power Source:** Operates on CO₂ or N₂O gas
(no electrical power required)

Cryo Carriers

- CR4150 Dual cylinder cryo carrier with three 30 inch hoses & one switchover valve (holds either 2 C or 2 E cylinders)

Cryo Filters

- CR4400 Microfiber filter retrofit
- CR4350 12 filters and 6 o-rings
- CR5350 Sintered metal filters (pack of three)

Gas Connectors

- CR4210 Handtight N₂O for C cylinder
- CR4220 Handtight CO₂ for C cylinder
- CR4280 T-Yoke N₂O for E cylinder
- CR4290 T-Yoke CO₂ for E cylinder
- CR4310 Handtight N₂O for C cylinder with filter
- CR4320 Handtight CO₂ for C cylinder with filter
- CR4380 T-Yoke N₂O for E cylinder with filter
- CR4390 T-Yoke CO₂ for E cylinder with filter

Cryo Accessories

- CR4160 25 foot exhaust hose
- CR4260 60 inch pressure hose
- CR4360 30 inch pressure hose
- CR4270 Protective probe tip cap