


LASER PROBES


While there are many similarities between laser probes these days, we have focused on several key performance areas to deliver a significant improvement to your endo-photocoagulation experience.



"Curve angle 90+ degrees, smooth actuation and perfect spot size...best probe on the market"
Surgeon comment

MAXReach LASER PROBE

Gauge	Product#	Description
25ga	VS0140.25	MAXReach
23ga	VS0140.23	MAXReach
25ga	VS0140.25E	MAXReach Extendable (tip advance)
23ga	VS0140.23E	MAXReach Extendable (tip advance)
25ga	VS0141.25	MAXReach Aligned (button and curve aligned)
23ga	VS0141.23	MAXReach Aligned (button and curve aligned)



"Smooth transition at taper, easy pass through cannula"
Surgeon comment

FLEX-TIP LASER PROBE

Gauge	Product#	Description
25ga	VS0130.25	25ga Flex-Tip
23ga	VS0130.23	23ga Flex-Tip
25ga	VS0135.25	25ga Illuminated Flex-Tip
23ga	VS0135.23	23ga Illuminated Flex-Tip



CURVED LASER PROBE

Gauge	Product#	Description
27ga	VS0120.27	27ga Curved
25ga	VS0120.25	25ga Curved
23ga	VS0120.23	23ga Curved



STRAIGHT LASER PROBE

Gauge	Product#	Description
27ga	VS0110.27	27ga Straight
25ga	VS0110.25	25ga Straight
23ga	VS0110.23	23ga Straight

All Vortex probes are supplied with the following premium features:

- Premium malleable fiber cladding eliminates "coiling" on the surgical field
- Hand polished fibers for premium finish and maximum output

Are vigorously tested for:

- Uniform spot size
- Efficiency of emission, defect and laser power output
- Smoothness of moving components, cannula access and laser connectivity

And are provided:

- Sterile, 5/box, SMA 905 COMPATIBLE

Compatible with:

Iridex, Alcon, Bausch & Lomb, Ellex, Zeiss, DORC (US model)

Trademarks are the property of their respective owners