


Retinal Cannulas



Product#	Description
VS0210.25	25ga Slotted Soft Tip
VS0210.23	23ga Slotted Soft Tip

ESTAFANOUS SLOTTED SOFT TIP

- "Slotted" tubing allows for easy transition thru valves.
- Particularly useful for detachments allowing complete drain of fluid while minimizing tissue incarceration.
- Rotate handle to allow for complete fluid removal.
- Sterile, 5/box



Product#	Description
VS0200.25	25ga Soft Tip Cannula
VS0200.23	23ga Soft Tip Cannula

SOFT TIP CANNULA

- Custom tubing provides
 - High flow
 - Rigid enough to go through valved cannulas
 - Soft enough for retinal surface
- Sterile, 5/box
- 1mm tubing length



Product#	Description
VS0280.25	25ga Universal "RETRACTABLE"
VS0280.23	23ga Universal "RETRACTABLE"

VERSATILE SOFT TIP "RETRACTABLE"

- Single use disposable
- .45in/11.5mm extended silicone tubing length
- Extendable silicone tubing allows for injection or drainage under the retina
- Soft, silicone tubing allows for atraumatic injection or drainage on the retina surface
- Sterile, 5/box



Product#	Description
VS0220.25	25ga Sub-Retinal Injection Cannula

SUB-RETINAL INJECTION CANNULA

- 25ga shaft with 41ga x 88mm long polyimide tubing
- Rigid polyimide allows easy perforation through retina
- Sterile, 5/box

