


E.D.D. - External Drainage and Depression

PRODUCT SPECIFICS:

Bulb depressor - .086 diameter

Needle Size: 28ga

Exposed Needle Length: 2.4mm



BENEFITS:

- ◆ Eliminates the need for a cut down
- ◆ Precise localization of 28ga drainage point
- ◆ Fixed drainage needle length for complete control
- ◆ Potentially reduces the risk of retinal incarceration
- ◆ Potentially reduces the risk of hemorrhages

Product Code: VS0290, Sterile, 5/bx

For more information, please contact Vortex Surgical at
636-778-4350 or info@vortexsurgical.com

