

Circling Bands

S 2969 (40 Style)
To be used with 1, 2, 4.

S 2987 (240 Style)
To be used with 1, 3, 4.

S 2950 (2950 Style)
To be used with 1, 3, 4.

Silicone Lace™

S 5020 (40 Style)
To be used with 1, 2, 4.

S 5025 (240 Style)
To be used with 1, 3, 4.

Silicone Strips

S 2970 (41 Style)
To be used with 5, 6, 7.

S 2971 (42 Style)
To be used with 5, 6, 7.

Grooved Strips

S 2992 (219 Style)

S 2968 (32 Style)

S 2998 (220 Style)

S 3013 (225 Style)

Silicone Sleeves

Round

S 3018 (70 Style)

S 3019 (270 Style)

S 3071 (72 Style)

Oval

S 3083 (3083 Style)

S 3084 (3084 Style)

Wong Buckling Sleeves™

S 3093 (3093 Style)

S 3094 (3094 Style)

Since 1959

LABTICIAN™

ophthalmics

RETINAL IMPLANTS

Silicone Sponges

Round

S 1982-2 (502 Style)

S 1982-3 (503 Style)

S 1982-4 (504 Style)

S 1982-5 (505 Style)

S 1982-7 (8270 Style)

Part Thickness

S 1984-2.5 (8425 Style)

S 1984-3.3 (515 Style)

S 1984-5 (510 Style)

S 1984-7.5 (511 Style)

Grooved

S 1983-2.3 (516G Style)

S 1983-2.5 (510G Style)

S 1983-3 (506G Style)

S 1983-3.5 (509G Style)

Oval

S 1981-4 (501 Style)

S 1981-5 (506 Style)

S 1981-7 (907 Style)

S 1981-7.5 (507 Style)

Scholda

S 1985-5 (8550 Style)

S 1985-7 (8570 Style)

Silicone Tires

Concave

S 2986 (277 Style)

S 2991 (279 Style)

Convex

S 3010 (286 Style)

S 3014 (287 Style)

Asymmetrical

S 2994 (275 Style)

S 2999 (276 Style)

S 2996 (280 Style)

Wide Strips

S 4050 (4050 Style)
To be used with 5, 6, 7.

S 4060 (4060 Style)

Accessories

Meridional Implants

S 2978 (103 Style)

Button

S 2976 (22 Style)

LINCOFF RULES FOR FINDING PRIMARY BREAK

98%

SUPERIOR TEMPORAL OR NASAL DETACHMENTS

In 98% of the cases, the primary break lies within 1-1/2 clock hours of the highest border

93%

TOTAL OR SUPERIOR DETACHMENTS THAT CROSS THE 12 O'CLOCK MERIDIAN

In 93% of the cases, the primary break is at 12 o'clock or in a triangle, the apex of which is at the ora serrata and the sides of which intersect the equator 1 hour on either side of 12 o'clock

95%

INFERIOR DETACHMENTS

In 93% of the cases, the higher side of the detachment indicates on which side of the 6 o'clock radian an inferior break lies. When an inferior detachment is bullous, the primary break lies above the horizontal meridian

Cross Reference Chart					
Style	Product	Implant Type	Style	Product	Implant Type
22	S 2976	Silicone Button	279	S 2991	Silicone Tire
32	S 2968	Grooved Strip	280	S 2996	Silicone Tire
40	S 2969	Circling Band	281	S 2997	Silicone Tire
40	S 5020	Silicone Lace	286	S 3010	Silicone Tire
41	S 2970	Silicone Strip	287	S 3014	Silicone Tire
42	S 2971	Silicone Strip	287 WG	S 3014 L	Silicone Tire
70	S 3018	Silicone Sleeve	501	S 1981-4	Silicone Sponge Oval
72	S 3071	Silicone Sleeve	502	S 1982-2	Silicone Sponge Round
103	S 2978	Meridional Implant	503	S 1982-3	Silicone Sponge Round
219	S 2992	Grooved Strip	504	S 1982-4	Silicone Sponge Round
220	S 2998	Grooved Strip	505	S 1982-5	Silicone Sponge Round
225	S 3013	Grooved Strip	505 G	S 1983-5	Silicone Sponge Grooved
240	S 2987	Circling Band	506	S 1981-5	Silicone Sponge Oval
240	S 5025	Silicone Lace	506 G	S 1983-3	Silicone Sponge Grooved
270	S 3019	Silicone Sleeve	507	S 1981-7.5	Silicone Sponge Oval
275	S 2994	Silicone Tire	509 G	S 1983-3.5	Silicone Sponge Grooved
276	S 2999	Silicone Tire			
277	S 2986	Silicone Tire			
510	S 1984-5	Silicone Sponge Part Thickness	907	S 1981-7	Silicone Sponge Oval
510 G	S 1983-2.3	Silicone Sponge Grooved	2950	S 2950	Circling Band
511	S 1984-7.5	Silicone Sponge Part Thickness	3083	S 3083	Silicone Sleeve Oval
515	S 1984-3.3	Silicone Sponge Part Thickness	3084	S 3084	Silicone Sleeve Oval
516 G	S 1983-2.3	Silicone Sponge Grooved	3093	S 3093	Silicone Sleeve Buckling
			3094	S 3094	Silicone Sleeve Buckling
			4050	S 4050	Silicone Strip Wide
			4060	S 4060	Silicone Strip Wide
			8270	S 1982-7	Silicone Sponge Round
			8425	S 1984-2.5	Silicone Sponge Part Thickness
			8550	S 1985-5	Silicone Sponge Scholda
			8570	S 1985-7	Silicone Sponge Scholda

Indicates products exclusive to Labtician.

Manufactured by Labtician Ophthalmics, Inc.
2140 Winston Park Drive, Unit 6, Oakville, ON L6H 5V5 Canada
Tel: (905) 829-0055 or toll free (800) 265-8391
Fax: (905) 829-0056 info@labtician.com